

PHCBI Multigas Incubator MCO-170M



161L; Dual IR and Stabilized Zirconia Sensor; 5-50C

Rating: Not Rated Yet

Price

Sales price \$11,950.00

[Ask a question about this product](#)

Manufacturer [Panasonic](#)

Description

The double-stack MCO-170M Multigas Incubator offers precision regulation of oxygen and other gases, while maintaining temperature uniformity, contamination control, and usability.

Per chamber:

External Dimensions (W x D x H): 620 x 730 x 905 mm

Internal Dimensions (W x D x H): 490 x 523 x 665 mm

Volume: 161 liters

Temperature Control Range & Fluctuation: AT +5 ~ +50, ± 0.1 °C

Temperature Uniformity ± 0.25 °C

CO₂ Control Range & Fluctuation: 0 ~ 20, ± 0.15 %

Incubators: PHCBI Multigas Incubator MCO-170M

O2 Control Range & Fluctuation: 1 -18 and 22 - 80, ± 0.2 %
Humidity Level & Fluctuation: 95, ± 5 %RH
Sterilisation Method: H2O2 Decontamination
Temperature Sensor: Thermistor
CO2 Sensor: Dual IR
O2 Sensor: Stabilized Zirconia Sensor
Inner Doors: 4 gastight - made of tempered glass
Shelves: 3 supplied as standard
Shelf Dimensions (W x D x H): 470 x 450 x 12 mm

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com