

## PHCBI Multigas Incubator MCO-170M



161L; Dual IR and Stabilized Zirconia Sensor; 5-50C

Rating: Not Rated Yet

**Price**

Sales price \$11,950.00

[Ask a question about this product](#)

Manufacturer [Panasonic](#)

### Description

The double-stack MCO-170M Multigas Incubator offers precision regulation of oxygen and other gases, while maintaining temperature uniformity, contamination control, and usability.

Per chamber:

External Dimensions (W x D x H): 620 x 730 x 905 mm

Internal Dimensions (W x D x H): 490 x 523 x 665 mm

Volume: 161 liters

Temperature Control Range & Fluctuation: AT +5 ~ +50,  $\pm 0.1$  °C

Temperature Uniformity  $\pm 0.25$  °C

## Incubators: PHCBI Multigas Incubator MCO-170M

---

CO2 Control Range & Fluctuation: 0 - 20,  $\pm 0.15$  %  
O2 Control Range & Fluctuation: 1 - 18 and 22 - 80,  $\pm 0.2$  %  
Humidity Level & Fluctuation: 95,  $\pm 5$  %RH  
Sterilisation Method: H2O2 Decontamination  
Temperature Sensor: Thermistor  
CO2 Sensor: Dual IR  
O2 Sensor: Stabilized Zirconia Sensor  
Inner Doors: 4 gastight - made of tempered glass  
Shelves: 3 supplied as standard  
Shelf Dimensions (W x D x H): 470 x 450 x 12 mm

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**