

## Precision 180 Series Water Bath

PRODUCT  
Image  
Unavailable

PRODUCT  
Image  
Unavailable

Ideal for applications in which temperature setpoints seldom change Variable analog temperature control Back-up high-limit safety thermostat Temperature displayed on 0 to 100°C glass thermometer Useful for a wide variety of application

Rating: Not Rated Yet

**Price**

[Ask a question about this product](#)

Manufacturer [Accuri](#)

### Description

- Ideal for applications in which temperature setpoints seldom change
- Variable analog temperature control
- Back-up high-limit safety thermostat
- Temperature displayed on 0 to 100°C glass thermometer
- Useful for a wide variety of applications

Exceptional price, performance and reliability; for most general laboratory, environmental, medical and biotechnology applications. Broad temperature range: ambient up to 99.9°C. Excellent temperature uniformity:  $\pm 0.2^{\circ}\text{C}$  at 37°C. Eight different sizes available. Easy to use and maintain. Temperature setpoint is easily adjusted/programmed from control panel.

Features over-temperature safety, fully insulated chamber with stainless steel diffuse shelf and gable cover. Bath temperature is displayed on a 0 to 100°C glass thermometer. One-piece stainless steel interior and epoxy-powder coated exterior make unit easy to clean.

1.5 Liter capacity (.4 gallons)