

Q Sonica Sonicator Sonic Dismembrator





Model Q55

Rating: Not Rated Yet

Price

Sales price \$1,250.00

[Ask a question about this product](#)

Description

Compact and cost-effective ultrasonic processor. Effective for standard cell disruption and many other small volume applications. 55W; 24VDC; 20kHz

Other specifications:

Dimensions: 8" W x 7.5" L x 5.75" H

Volume: Up to 50 ml

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com