

Rainin Pos-D Positive Displacement Pipettors -- New in Box!



Most sizes are available.

Rating: Not Rated Yet

Price

Sales price \$250.00

[Ask a question about this product](#)

Description

New in box!

Manual pipettor.

Positive displacement pipette using piston-in-capillary tips for precise, contamination-free transfer of viscous or volatile liquids.

Most sizes available.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.

South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com