

Retsch MM2 Roller Mill Shaker



Variable speed mixer/shaker.

Rating: Not Rated Yet

Price

Sales price \$1,500.00

[Ask a question about this product](#)

Manufacturer [Retsch](#)

Description

- Fully functional with Variable speed.
- 2 holders 1 1/4 Diameter (Opens to 2 1/2? Adjusts down to 2 1/8?)
- Perform radial oscillations in a horizontal position.
- intensive mixing of the sample
- Application Examples: alloys, animal feed, bones, ceramics, cereal grains, chemical products, coal, coke, drugs, glass, hair, minerals, oil seeds, ores, paper, plant materials, plastics, sewage sludge, soils, tablets, textiles, tissue, waste samples, wood, wool, ...
- Size reduction, mixing, homogenization, cell disruption
- Field of application: agriculture, biology, chemistry / plastics, construction materials, engineering / electronics, environment, food, glass / ceramics, medicine / pharmaceuticals, mineralogy / metallurgy
- Feed material: hard, medium-hard, soft, brittle, elastic, fibrous, Size reduction principle impact, friction