

REVCO Elite II



Revco Elite II double-stack, water jacketed CO2 incubator w/IR sensor

Rating: Not Rated Yet

Price

Sales price \$6,500.00

[Ask a question about this product](#)

Manufacturer [REVCO](#)

Description

!!! SOLD !!!

This Revco Elite II double-stack, water jacketed CO2 incubator, model RCO5000D-5-ABC, 2 x 5.8 cu/ft chambers, infrared (IR) sensor, microprocessor IntrLogic II controlled, user-selectable RH display with 1% resolution, water jacket temp indicator, 1.5" (38 mm) porthole (in rear wall), is very clean and lightly used.

Please call us at (650) 871-5707 to check availability of this product!