

## Sanyo MCO-17AIC



Sanyo MCO-17AIC double-stack CO2 incubator w/IR sensor

Rating: Not Rated Yet

**Price**

Sales price \$4,200.00

[Ask a question about this product](#)

Manufacturer [Sanyo](#)

Description

**!!!! - RECENTLY SOLD - !!!!**

This Sanyo MCO-17AIC double-stack CO2 incubator unites stainless steel (corrosion resistance) with copper (antibacterial) and provides you with a sterile interior made of electro-polished copper alloy stainless. 2 x 5.8 cu ft chambers, direct heat and air jacket system, 5 to 50 C, 0 to 20% CO2