

Sartorius Entris II 6202i Precision Balance



Maximum capacity: 6200 g

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Sartorius](#)

Description

Sartorius Entris II 6202i precision balance has internal isoCAL adjustment. Max 6200 g, d = 10 mg, and e = 100 mg, with a CE M21 mark. It has a stainless-steel pan, a hybrid LED/touch display, and the orange tare key typical of the Entris II line.

Key specifications

Maximum capacity: 6200 g

Readability (d): 10 mg (0.01 g); verification interval (e): 100 mg

Typical stabilization time: ? 0.9 s (typical repeatability about 5 mg at 5% load; linearity ± 20 mg)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com