

Sartorius Practum 612-1S Precision Balance



Maximum capacity: 610 g

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Sartorius](#)

Description

The Sartorius Practum 612-1S (often written PRACTUM612-1S) is an entry-level precision balance for routine lab weighing. It has a color touchscreen with swipe operation, external calibration, a USB interface, overload protection, and built-in applications such as weighing, counting, percent weighing, and check weighing.

Key specifications

Maximum capacity: 610 g

Readability: 0.01 g (10 mg)

Typical stabilization time: 1.5 s (repeatability 0.01 g; linearity deviation 0.03 g)

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com