

Sartorius SECURA 612-1S Top-Loading Balance



Capacity x Readability: 610 x 0.01 g

Rating: Not Rated Yet

Price

Sales price \$975.00

[Ask a question about this product](#)

Manufacturer [Sartorius](#)

Description

A high-quality, reliable top-loading balance.

Technical specifications

Repeatability: 10mg

Linearity: 20mg

Minimum sample weight according to USP, Typical: 12g

Temperature limits: or verified models with activated isoCAL +10°C | +30°C

Typical stabilization times (seconds): 1.5g

Weighing pan size, diameter: 180mm

Dimensions (DxWxH): 360 x 216 x 95 mm - 14.1 x 8.5 x 3.75 inch

Balances: Sartorius SECURA 612-1S Top-Loading Balance

Levelling: Intelligent, optoelectronic level sensor with alarm function and interactive user guidance for reliable and safe leveling

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com