

Sartorius YDP03-0CE Printer



Sartorius YDP03-0CE data printer with built-in real-time clock with statistics programs designed for Sartorius balances equipped with an RS-232 interface.

Rating: Not Rated Yet

Price

Sales price \$450.00

[Ask a question about this product](#)

Manufacturer [Sartorius](#)

Description

Text and data are printed with up to 24 characters per line and the unit is equipped with a real-time clock. High contrast LCD that shows the date or time as well as the current transaction number. Integrated rechargeable batteries with battery level indicator. Runs on the balance's power, no separate power supply necessary. Easy-to-change paper roll and ink ribbon cassette.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com