

## Savant Refrigerated VaporTrap RVT400



Effectively condenses and traps a wide variety of solvent vapors that can contaminate and degrade vacuum pump oil. Easy-to-use glass condensation flask (SAV-GCF400) allows for the collection of liquid into a self contained jar for easy draining. Rapid cool

Rating: Not Rated Yet

**Price**

Sales price \$800.00

[Ask a question about this product](#)

Manufacturer [Savant](#)

### Description

Effectively condenses and traps a wide variety of solvent vapors that can contaminate and degrade vacuum pump oil. Easy-to-use glass condensation flask (SAV-GCF400) allows for the collection of liquid into a self contained jar for easy draining. Rapid cooling time minimizes waiting period.

Temperature -50c

Capacity 4 liters

Dimensions WxHxD 14x12x24 in

Weight, lbs 55