

Savant UVS400



SOLD The Savant UVS400 Universal Vacuum System combines an oil-free vacuum source and refrigerated vapor trap in one.

Rating: Not Rated Yet

Price

Sales price \$2,250.00

[Ask a question about this product](#)

Manufacturer [Savant](#)

Description

SOLD

Simple automated operation. Includes Glass Condensation Flask for rapid emptying of the solvent without the need to wait until the solvent defrosts.

Small footprint helps maximize valuable bench space.

Designed for use with concentrators, gel dryers and any operation requiring vacuum and solvent recovery.

Universal Vacuum Systems (UVS) simplify sample preparation and setup by combining the vacuum source and the refrigerated vapor trap into a single instrument. There is less tubing to plumb, reduced maintenance and best of all, it can be used for almost any application requiring vacuum and solvent recovery.