

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Seward Stomacher 400



Improved organism recovery with the patented curved shape of the Circulator paddles over older square paddle designs No risk of cross contamination between samples Durable and reliable Control panel with programmable settings or 'hands-free'

Rating: Not Rated Yet

Price

Sales price \$1,900.00

[Ask a question about this product](#)

Manufacturer [Seward](#)

Description

WHY the Stomacher® 400 Circulator

- Improved organism recovery with the patented curved shape of the Circulator paddles over older square paddle designs
- No risk of cross contamination between samples
- Durable and reliable
- Control panel with programmable settings or 'hands-free' auto-run
- Countdown of processing time
- Removable door for easy cleaning of paddle chamber

Meeting the food safety and quality demands of today's microbiologists.

Whatever method you use to enumerate and identify bacteria, yeasts and moulds in your laboratory, thorough reproducible homogenisation of the sample is the first and most critical step.

Numerous studies have demonstrated the Circulator enhancement as a positive benefit to pathogen testing procedures in traditional microbiology and PCR based techniques.

As the sample is crushed during operation, the patented curved shape of the Circulator paddles create a washing effect driving more organisms in suspension thereby providing a better homogenised sample for analysis.

The Stomacher® 400 Bag Range

Caters for both the [Circulator bag](#) (round-bottom) designed bags, and the [Classic](#) (square-bottomed) for the older square paddle designed

Stomachers®.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

[Biotech Equipment Sales, Inc.](#)
[226 Miller Ave.](#)
[South San Francisco, CA 94080](#)
[650-871-5707 office](#)
[650-276-7487 fax](#)
sales@biotechequipmentsales.com