

Sorvall Biofuge Fresco Refrigerated Centrifuge



The compact refrigerated Biofuge Fresco microlitre centrifuge creates space in the laboratory. Even at maximum speed it operates extremely quietly, Short acceleration and braking times ensure rapid sample processing.

Rating: Not Rated Yet

Price

Sales price \$2,000.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

!!!SOLD!!!

The compact refrigerated Biofuge Fresco microlitre centrifuge creates space in the laboratory. Even at maximum speed it operates extremely quietly, Short acceleration and braking times ensure rapid sample processing.

The brushless induction motor of the Fresco is extremely efficient. In 15 seconds, it accelerates the rotor to its maximum RCF of 16,060xg, which is sufficient for all standard applications. A "quick run" key is used for short processing runs.