

Sorvall F-16/250 Fixed-Angle Rotor



Sorvall F-16/250 Fixed-Angle Rotor for use in Sorvall Evolution RC centrifuge

Rating: Not Rated Yet

Price

Sales price \$1,500.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

Rotor F-16/250*

Capacity(piece x ml) - 6 x 250

Max Speed (rpm) - 14,500 in Evolution RC centrifuge

Max.RCF (x g) - 31,940 in Evolution RC centrifuge

Accel/Decel (min:sec) - 2:30/1:45

Tube Size (range, ral) - 3-250

Catalog Number - 11966

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com

