

Sorvall RC 3C Plus



Sorvall RC 3C Plus digital refrigerated floor model centrifuge w/six by 1 liter rotor

Rating: Not Rated Yet

Price

Sales price \$8,500.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

This Sorvall RC 3C Plus General Purpose digital refrigerated floor model centrifuge w/the H6000A six by 1 liter rotor is ideal for clinical labs, blood banks, or industrial settings. It is microprocessor controlled and has programming capabilities along with other automated features. The temp range is 0 to 30C and it can get to 5000 rpm w/this rotor.

The centrifuge is pristine. The rotor shield shows some corrosion, as well as some spots on a few of the buckets, but the rotor body is very clean. Current model, manufactured in 1997.