

## Sorvall RC 5B



Sorvall RC 5B high-speed refrigerated floor model centrifuge

Rating: Not Rated Yet

**Price**

Sales price \$4,750.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

**!!!! - RECENTLY SOLD - !!!!**

This Sorvall RC 5B high-speed refrigerated floor model centrifuge gets to 20,000 rpm, can accept a variety of rotors for different sized tubes, temp range from -20 to 40C, max volume is 3 liters (6 x 500 ml)