

Sorvall SA-600



Sorvall SA-600 Fixed-Angle Rotor, 12 Places, 50mL each

Rating: Not Rated Yet

Price

Sales price \$1,300.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

SORVALL SA-600 Fixed-Angle Rotor and lid in good condition. For use with Sorvall Evolution RC and RC-5 series centrifuges.

Specifications:

12 places, 50 mL ea.

Diameter: 27cm (10-3/4 in)

Mass (weight): 10.2 kg (22.5 lbs)

Angle: 34 degrees

Number of Places: 12

Maximum Angular Velocity (in rev/min): 16,500/min

Maximum Relative Centrifugal Force (RCF): 39,400

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080

650-871-5707 office
650-276-7487 fax
sales@biotequipmentsales.com