

Sorvall SA 600 Fixed-Angle Rotor



For use with Thermo Scientific™ Sorvall™ RC 6 Plus and Evolution RC superspeed Centrifuges.

Rating: Not Rated Yet

Price

Sales price \$1,950.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

Processes samples efficiently.

Material: Aluminum

Capacity: 12 x 50 mL

Maximum Speed 17,000 rpm

Autoclavable 121 C

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com