

Sorvall SLA 3000 Super Lite Rotor



Fixed-angle rotor, max speed 12,000 RPM.

Rating: Not Rated Yet

Price

Sales price \$3,250.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description

The Sorvall SLA 3000 Super Lite Rotor has a capacity of 6 x 500 ml, along with 6 smaller tube holders. Max RCF of 24,336 x g and autoclavable at 121°C. The SLA 3000 is compatible with Thermo Scientific Sorvall Evolution RC series and RC-6 Plus series centrifuges (including RC-5, RC-5B, RC-5B Plus, RC-5C, RC-5C Plus, and RC-26 Plus). K factor: 2364.

[Biotech Equipment Sales, Inc.](#)
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com