

Sorvall SS-34 Fixed Angle Rotor



Multi-liter cell culture harvesting or the collection of cell-free growth media containing excreted products such as proteins, hormones and vitamins.

Rating: Not Rated Yet

Price

Sales price \$1,500.00

[Ask a question about this product](#)

Manufacturer [Sorvall](#)

Description Specifications

Type	Fixed-Angle
Material	Aluminum
Capacity	8 x 50mL
Maximum Speed	20,500rpm
Tube Dimension (O.D.xL)	29 x 108mm
Angle	34°
K Factor	714
Net Weight	6.8kg (15.0 lb.)
RCF-max	50,228xg
RCF-min	15,350xg
Radius-max	107.0mm
Radius-min	32.7mm

Biotech Equipment Sales, Inc.

226 Miller Ave.

South San Francisco, CA 94080

650-871-5707 office

650-276-7487 fax

sales@biotechequipmentsales.com