

## SP Scientific GeneVac EZ-2.3 Plus Solvent Evaporating System

2000 rpm maximum speed, 500xg maximum g-force,

Rating: Not Rated Yet

**Price**

Sales price \$11,500.00

[Ask a question about this product](#)

### Description

This version is HCl compatible, Recommended for difficult solvents with higher boiling points up to 165°C, suitable for water and HPLC fractions, features 2000 rpm maximum speed, 500xg maximum g-force, 0-1200 mbar vacuum pressure display, automatic/2 mbar - atmosphere pressure control, 2 mbar system ultimate pressure, ambient to 60°C chamber temperature range, -50°C minimum condenser temperature, solvent compatibility includes alcohols, DCM/methylene chloride, ethyl acetate, water, TFA, DMF, HCl, and HNO<sub>3</sub>, a dual-place swinging bucket rotor is included, a diaphragm vacuum pump is integrated inside of the system, 120V.

\*Two sets of sample holders are included; for 16mm Ø tubes and 20mL scintillation vials. Other sample holders may be available, please inquire for details.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)