

Stemcell Technologies EasySep Magnet



Designed to hold a standard 12 x 75 mm (5 mL) polystyrene tube.

Rating: Not Rated Yet

Price

Sales price \$500.00

[Ask a question about this product](#)

Description

Magnet for column-free immunomagnetic separation

Easily and efficiently perform magnetic cell separation procedures by using the EasySep Magnet with EasySep cell isolation reagents. The EasySep Magnet generates a high-gradient magnetic field in the interior cavity that is strong enough to separate cells labeled with EasySep magnetic particles without the use of columns. This magnet is designed to hold a standard 12 x 75 mm (5 mL) polystyrene tube.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com