

Stemi 2000C Stereo Trinocular Microscope with Illuminator



The Zeiss Stemi 2000-C stereo microscope offers a wide variety of standard features and options that will meet the most rigorous application demands

Rating: Not Rated Yet

Price

Sales price \$1,850.00

[Ask a question about this product](#)

Manufacturer [Stemi](#)

Description

The sharp, distortion-free 3D images produced by the Greenough optical system, along with a 7.7X to 1X click stop or continuous zoom range, sets the Stemi 2000-C Zeiss Microscope apart from other stereo microscopes. The included 10X eyepieces provide a 6.5X to 50X magnification range yielding a field of view from 35.4mm down to 4.6mm. By adding supplemental lenses and 25X eyepieces, magnifications of 1.95X to 250X can be achieved for this Zeiss Microscope.

- 7.7X to 1X Zoom Range
- Video Adaptable
- Long Working Distance
- Easy Boom Stand Adaptation
- Magnifications of 1.95X to 250X
- Field of View from 118mm to 1.0mm