

## Stemcell Technologies The Big Easy EasySep Magnet



Designed to hold a standard 14 mL (17 x 95 mm) polystyrene tube.

Rating: Not Rated Yet

**Price**

Sales price \$650.00

[Ask a question about this product](#)

### Description

Magnet for column-free immunomagnetic separation

Easily and efficiently perform magnetic cell separation procedures using "The Big Easy" EasySep Magnet with EasySep cell isolation reagents. This magnet is also used as a component of RoboSep instruments when processing large numbers of cells. "The Big Easy" EasySep Magnet generates a high-gradient magnetic field in the interior cavity that is strong enough to separate cells labeled with EasySep magnetic particles without the use of columns. "The Big Easy" EasySep Magnet is designed to hold a standard 14 mL (17 x 95 mm) polystyrene tube.

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)