

Thermo Biosafety Cabinet 1345



4 ft.; with motorized lift stand.

Rating: Not Rated Yet

Price

Sales price \$6,250.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

Class II Type A2 (70% air recirculated and 30% exhausted from a common plenum to the room). Comes with motorized lift stand. 120V.

MULTIPLE UNITS AVAILABLE.

Comes with our 90-day warranty for all parts, electrical and mechanical, except the HEPA filters, which are guaranteed to pass initial certification or BES will pay for 50% of replacement cost.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com