

## Thermo Compact Digital Microplate Shaker, 4-place



4-place

Rating: Not Rated Yet

**Price**

Sales price \$525.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

Catalogue # 88880023

Model: Compact Digital Microplate Shaker

Capable of shaking 4 plates in one run.

Max load: 8.8 lbs.

Electrical: 100-240V 50/60Hz

Features:

PID control ensures consistent shaking results

Digital display for easy viewing

Compact design is ideal for placing in incubators (non-CO<sub>2</sub>), chambers, and refrigerators (up to 60°C)

Timer Range: Continuous or 1min to 99hrs 59min

Automatic Overload Shutdown: 4.4kg (9.7lb.)

Small 3mm orbital motion, ideal mixing small volumes

Tube rack accessories make this unit extremely versatile

Includes 4 place microplate

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**