


## Thermo Flex 5000 Recirculating Chiller



208-230V

Rating: Not Rated Yet

**Price**

Sales price \$6,250.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

Functional and efficient, with an internal bath capacity of 7.2 liters, and a flow rate of 4.1 GPM at 60 Psi. This robust recirculating chiller also includes an LCD digital display that presents current setpoint temperature settings and actual solution temperature. 208-230V

Temp range: 5 to 40C

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)