

Thermo Forma 906



Thermo Forma 906 23 cu/ft upright -80 ultralow freezer. Deluxe series. 208V.

Rating: Not Rated Yet

Price

Sales price \$5,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Forma](#)

Description

!!!SOLD!!!

This Thermo Forma 906 23 cu/ft upright -80 ultralow freezer. Deluxe series. -5-C to -86C. The centralized information center, which includes a microprocessor control and monitoring system and an optional chart recorder or datalogger, ensures that all controls and displays are easy-to-reach and read. The 900 series is designed for universities, research facilities and industrial institutions where easily adjustable shelves, easy-to-access controls and convenient at-a-glance monitoring are preferred. 208 V, 12A