

## Thermo Forma Series 3 Doublestack Water-Jacketed CO2 Incubator





6.5 cubic foot chamber (184L)

Rating: Not Rated Yet

**Price**

Sales price \$8,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Forma](#)

Description

Model 4110

The Thermo Forma 4110 is a water-jacketed unit with a 6.5 cubic foot chamber, polished stainless steel interior. Unique triple wall construction provides outstanding temperature uniformity and stability. The heated dual pane inner glass door is designed to minimize condensation. An in-chamber fan gently distributes air throughout the chamber, designed to ensure homogeneous conditions, while the in-chamber HEPA air filtration continuously filters the entire chamber volume every 60 seconds. The stainless steel shelves and supports can be removed without tools for easy cleaning or adjustment.

Each Thermo Forma Series 3 4110 water-jacketed CO2 incubator has:

6.5 cubic foot chamber (184L)

Temperature range: ambient +5°C to 50°C

T/C CO2 sensor

Thermo iCAN touch screen sensor

Specifications

Capacity (English): 6.5 cu. ft.

Certifications/Compliance: UL, cUL, CE

Data Outputs: USB, and optional 4-20mA

Description: Single 184L incubator

Capacity (Metric): 184 L

Chamber Material: Polished Stainless Steel

CO2 Concentration Range: 0 to 20%

CO2 Sensor Technology: TC Sensor

Lighting: LED

Material: Stainless Steel

Oxygen Control: None

Relative Humidity: <90% at 37°C

Temperature Range (Metric); Ambient +5° to 50°C

Type: CO2 Incubator

Dimensions (D x W x H) Interior: 20.0 x 21.3 x 26.8 in. (50.8 x 54.1 x 68.1 cm)

Dimensions (L x W x H): 25 x 26 x 39.5 in. (63.5 x 66 x 100.3 cm)

Weight (English): 365 lb.

Electrical Requirements: 115 V, 50/60 Hz

Unit Size: Each

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)