

## Thermo Helios Aquamate 2000E



Aquamate



The Thermo Helios Aquamate 9423 AQA 2000E UV-Vis spectrometer

Rating: Not Rated Yet

**Price**

Sales price \$2,500.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

The **Thermo Helios Aquamate** 9423 AQA 2000E UV-Vis spectrometer with the extended 190-1100nm range of the UV-Visible version provides additional flexibility to analyse detergents and screen for organics. AquaMate provides optical performance previously unheard of in a water analysis spectrophotometer. The 2nm bandwidth optical system means AquaMate has the analyte selectivity and sensitivity to handle a wide range of environmental applications other than just water analysis using test kits. For example, AquaMate is probably the only water analysis spectrophotometer capable of retrieving useful information from the sharp peaks found in chlorophyll spectra or from organic contaminants that absorb in the ultraviolet region of the spectrum. Superb linearity, stability and signal-to-noise ratio also means you can trust AquaMate to take stable measurements from highly absorbing or light scattering samples. See [here](#) for complete details.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)