

Thermo HERAcCell Vios 160i Double-Stack CO2 Incubator



On roller base. From 2021!

Rating: Not Rated Yet

Price

Sales price \$9,650.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

Features Direct Heat Chamber (165L electropolished stainless steel chamber with compact footprint, easily stackable, plus adjustable, perforated shelving and reversible exterior door for added flexibility); THRIVE Active Airflow; ISO 5 HEPA Filtration; Protected Humidity Reservoir; In situ sensors and probes; Dual temperature probes for over temperature protection with operational back up; Steri-Run, high temperature deflection cycle; Intuitive, interactive iCAN™ Touchscreen Control.

Specifications

Incubators: Thermo HERAcCell Vios 160i Double-Stack CO2 Incubator

Power: 120V, 50/60Hz
Capacity: 5.8 cu. ft. (165 L/chamber)
Certifications/Compliance: CE, CSA
Data Outputs: USB, and optional 4-20mA
Chamber Material: Electropolished Stainless Steel
Dimensions: 18.5" W x 22.7" D x 23.9" H (47 x 57.6 x 60.7 cm) per chamber
Dimensions: 25.1" W x 34.6" D x 35.4" H (63.7 x 88 x 90 cm) per chamber
Temperature Range: Ambient +3°C to 55°C
Relative Humidity: <93% at 37°C
Material: Stainless Steel
Humidity Source: Covered integrated water reservoir
Weight: 183 lbs. (83 kg) per chamber
CO2 Sensor Technology: TC Sensor
Lighting: LED
Type: CO2 Incubator
Oxygen Control: None
CO2 Concentration Range: 1 to 20%
Catalog number: 50144906

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com