

## Thermo IEC 8850 8x50mL Rotor Multi/Multi RF



Thermo IEC 8850 8 x 50mL with rotor lid for Thermo IEC Multi/Multi RF Centrifuge

Rating: Not Rated Yet

**Price**

Sales price \$1,400.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

Thermo IEC 8850 Fixed-Angle Rotor - for use in IEC Multi and Multi RF Centrifuges. Used, like new condition.

### Specifications -

Capacity - 8 x 50ml (Falcon/Corning)

Max. Speed - 15100 rpm

rpm at 4C - 14700 rpm

Max Force - 25250 xg

rcf at 4C - 23900 xg

Radius - 9.9 cm

Accel Time - 95 sec

Decel Time - 60 sec

Advantages - Aerosol Containment Adapters for 50mL Round bottom and 15mL conical.

Biotech Equipment Sales, Inc.  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax

[sales@biotequipmentsales.com](mailto:sales@biotequipmentsales.com)