

## Thermo Labsystems Multidrop DW



Dispenses liquids into 96 standard and DW plates across 20 to 995 $\mu$ L volume range.

Rating: Not Rated Yet

**Price**

Sales price \$3,200.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

With quick plate and volume selection combined with easy maintenance, the Multidrop DW makes assay setup quick and easy. Precise dispensing ensures high-quality assay results and reproducible data; wide volume range provides flexibility and increased speed in assays requiring larger volumes.

- Dispenses liquids into 96 standard and DW plates across 20 to 995 $\mu$ L volume range.
- Convenient, autoclavable dispensing cassette
- Uses an 8-channel detachable and autoclavable dispensing cassette to ensure sterile conditions and to avoid cross contamination
- All reagent lines can be backflushed into the reagent bottle, minimizing the loss of expensive reagents

[Biotech Equipment Sales, Inc.](#)  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)