

Thermo Lindberg Blue M BF51732PBFMC-1 Box Furnace



0.6 cu.ft. capacity; 100-1200°C temperature range

Rating: Not Rated Yet

Price

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

- 0.6 cu.ft. capacity
- 100-1200°C temperature range
- features 25 program 500 segment control
- over-temperature protection
- flow meter
- vertical lift door
- 208/240V, 50/60Hz, 16-19A, 4500W

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-

quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com