

## Thermo MaxQ 4000 Incubator Shaker - High Temp



Goes to 80C.

Rating: Not Rated Yet

**Price**

Sales price \$4,250.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

The MaxQ 4000 Incubator Shaker SHKE4000 offers reliable orbital shaking partnered with digitally controlled temperatures with an adaptable interior. This benchtop orbital shaker roomy stainless steel perforated 18 x 18 platform can be customized with various clamps to suit needs and can handle six 2-liter flasks at the same time, supporting a max load of 50 lbs.

The variable speed of the Maxq4000 can reach up to 400 RPM powered by a triple eccentric drive DC brushless motor, agitating items in a 19 mm / 0.75" orbit. The easy-to-use SHKE4000 incubated shaker has a chamber temperature range of 10°C above ambient to 60°C, with over-temperature safety, a clear scratch-resistant lid for easy sample viewing, and a timer function, up to 60 minutes or in continuous operation all controlled from the upfront dials. For safety, these Thermo incubator shakers will power down when the lid is opened.

Model number: SHKE4000-5

## Incubator Shakers: Thermo MaxQ 4000 Incubator Shaker - High Temp

---

Specifications:

Speed Range: 15 to 500rpm

Timer: Continuous, timed .01-999hours or 0.1- 999 minutes

Temperature Range: Ambient plus 10°C to 80°C

Temperature Accuracy: +/- 0.1°C

Maximum Load Capacity; 50lbs, Six 2L flasks

Orbit: 1.9cm (0.75 in.)

Platform Dimensions: 18 x 18 "

Dimensions (L x W x H): 23x33x32"

Power: 120V / 60Hz / 5amp

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)