

## Thermo MaxQ 4000 Model 4342 Refrigerating Benchtop Incubator Shaker



Speed Range: 15 to 500 RPM,  $\pm 1$  RPM; Temp Range:

Rating: Not Rated Yet

**Price**

Sales price \$4,500.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

### Description

Model 4342

Thermo MaxQ 4000 Model 4342 Refrigerating Benchtop Incubator Shaker

Versatile, incubated and refrigerated, these shakers are ideal for cell culture, solubility studies, extraction procedures, protein expression and many other applications. They are spacious chambers in a footprint designed for the benchtop.

### Specifications:

Temperature Range: 15°C below ambient to 60°C

Outside Dimensions: 34 x 22.5 x 22" (86.4 x 57.2 x 55.9 cm)

Pan Size: 18 x 18" (45.7 x 45.7 cm)

Power: 120V / A / W

Type: Digital Refrigerated Shaker

Accuracy:  $\pm 0.1^\circ\text{C}$  at 37°C in flask

Orbit Diameter: 0.75" / 1.9 cm

Maximum Load: 50 lb / 22.7 kg

Speed Range: 15 to 500 RPM,  $\pm 1$  RPM

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**