

## Thermo Orion Advanced Conductivity Meter, Model 162A



Used for a wide variety of applications.

Rating: Not Rated Yet

### Price

Sales price \$150.00

[Ask a question about this product](#)

Manufacturer [Thermo](#)

Description

Can be used for a wide variety of applications: Water, acids, bases and other samples can be easily analyzed for conductivity, salinity, resistivity and total dissolved solids. Features microprocessor design, which automates complicated and time consuming calibration and measurement procedures for a wide variety of applications. The Orion DirectCal calibration technique allows calibration by direct input of actual calibration standard values, eliminating iterative adjustments of cell constants.

Features:

Large, bright LED displays Conductivity, Resistivity, TDS, Salinity results along with temperature and current cell constant.

Automatically compensates for temperature changes based on the temperature coefficient and reference temperature.

Data log stores up to 200 data points; including measurement value, temperature and date/time of logging.

Autosampler mode.

Calibration selection (auto, manual, direct or multiple point)

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)