


Thermo Precision Shaker/Bath



Thermo Scientific Precision Reciprocating Shaker Baths are designed for life science and QA/QC applications. Push-button input of desired temperature and shaking speed help to make the baths easy to use, while meeting the stringent standards for laborato

Rating: Not Rated Yet

Price

Sales price \$1,250.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

!!!! - RECENTLY SOLD - !!!!

- Microprocessor control with digital LED display and shaker control
- Set temperature (in increments of 0.1°C) and oscillation speed via control-panel touchpad for entry of temperature and shaking speed
- Shaking frequency: 30 to 200 oscillations/min.
- Stroke length adjustable at 0.5, 1.0 or 1.5 in. (13, 25 or 38mm)
- Microprocessor-controlled proportional integral temperature control and solid-state sensing probe provide excellent temperature sensitivity
- Overtemperature protection capability
- Automatic motor shutdown capability
- With adjustable high-limit thermostat and indicator light
- Rugged stainless-steel interior with easy to clean epoxy powder-coated exterior