


Thermo Rocker Incubator 4637



Benchtop; 14" x 14" platform; temp range between 5C and 70C.

Rating: Not Rated Yet

Price

Sales price \$1,175.00

[Ask a question about this product](#)

Description

Consistent heating and rocking. Easy to use with digital temperature control and dial-in speed control.

- Adjustable angle
- Rocking speed of 0 to 100 cycles per minute.
- 14" x 14" platform (as pictured) accommodates a variety of containers.
- Temperature range of between 5°C above ambient to 70°C monitored by electronic proportional temperature control.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com