

## Thermo Savant RVT400 Refrigerating Vapor Trap



-55C operating temp; 5L trap capacity

Rating: Not Rated Yet

**Price**

Sales price \$1,250.00

[Ask a question about this product](#)

Manufacturer [Thermo Savant](#)

### Description

RVT400 is a compact benchtop refrigerated vapor trap for solvent vapor condensation and collection, used with Thermo Savant SpeedVac concentrators and gel dryers, -55°C operating temperature, 4-Liter trap capacity

115V; 4A; 50/60 Hz

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**

650-871-5707 office

650-276-7487 fax

[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)