

Thermo Savant SC250EXP



Thermo Scientific Savant SC250EXP large-capacity benchtop SpeedVac Concentrator has everything you need to effectively concentrate high boiling point solvents and strong acids and bases.

Rating: Not Rated Yet

Price

Sales price \$5,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Savant](#)

Description

Thermo Scientific Savant SC250EXP large-capacity benchtop SpeedVac concentrator.

- 4 IR heat lamps for drying high boiling point solvents such as DMSO and DMF
- Dual vacuum control of RAMP and LEVEL eliminates samples from bumping and/or freezing, provides optimisation for various sample types
- FFKM[®] and PTFE coated materials along the entire vapor path for maximum chemical resistance and durability
- Dual valve mechanism for improved durability and exhausting
- Comes complete with enlarged vacuum port, extra wide tubing and new extra large mouth glass condensation flask
- Heated glass lid reduces drying times and prevents condensation
- Digital display of parameters (time, temperature & vacuum)
- Standard glass safety lids
- Solvents and Solvent Combinations - SC250EXP P1: Acetic acid, Acetone, Ammonium hydroxide, Formic acid, Hydrochloric acid, TFA, Methylene chloride and water, Drug metabolites in SPE/urine, DCM and heptane, Methylene chloride and MEOH, THF/Pyridine (95:5), DMSO and methanol
- Solvents and Solvent Combinations - SC250EXP P2: Chloroform, DMSO, DMF, Ethyl acetate, Hexane, Toluene, Methanol and chloroform

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com