

## Thermo Scientific Pyro-Multi-Magnestir



These Thermo Scientific Barnstead/Lab-Line Pyro-Multi-Magnestir Hot Plate Stirrers are for heating and stirring up to six flat-bottom vessels simultaneously. A solid-state speed control allows specific regulation of stirring to suit requirements. Hot plat

Rating: Not Rated Yet

**Price**

Sales price \$800.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

### Description

Thermo Scientific Pyro-Multi-Magnestir -

These **Thermo Scientific Barnstead/Lab-Line Pyro-Multi-Magnestir Hot Plate Stirrers** are for heating and stirring up to six flat-bottom vessels simultaneously. A solid-state speed control allows specific regulation of stirring to suit requirements. Hot plate heats quickly and uniformly to 320°C (608°F), with stepless thermostat. Top plate is non-warping, cast aluminum. Unit accommodates up to six 1000mL beakers at one time. Stirring range is up to 1200 rpm for liquids with viscosities from 0 to 100% glycerol. Independent controls enable the operator to stir, heat, or stir and heat liquids. Separate on-off switches allow stopping and restarting either heater or stirrer at any time without changing control settings. The power switches have individual pilot lights. The stirrer unit is powered with six permanently lubricated and fan-cooled motors. Housing is insulated to prevent table-top scorching.

Biotech Equipment Sales, Inc., based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
**226 Miller Ave.**  
**South San Francisco, CA 94080**  
**650-871-5707 office**  
**650-276-7487 fax**  
**[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)**