

## Thermo Scientific CryoMed Controlled Rate Freezer 7454



-180°C to +50°C; 1.7 cu ft.

Rating: Not Rated Yet

**Price**

Sales price \$11,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Model 7454

With microprocessor. Customizable. Features alarms and monitors.

(Requires LN2 tank that is not included.)

Cabinet Material: Stainless steel with powder-coat paint finish

Programmability: Yes

Sample Loading: Front

Temperature Range: -180°C to +50°C

Capacity: 1.7 cu ft.

Power: 120 V / 60 Hz /12 A

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)