

Thermo Scientifici CryoPlus1 Model 7400



The ideal combination of liquid nitrogen, storage reliability and microprocessor technology.

Rating: Not Rated Yet

Price

Sales price \$4,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

The CryoPlus Series LN2 Storage system provides precise and accurate control over all parameters without complicated programming.

- Diameter (English): 16 in.
- Capacity (Metric): 90L
- Dimensions (L x W x H) Exterior: 26 x 21.5 x 41 in. (66 x 54.6 x 104cm)
- Type: Auto-Fill LN2 Cryo
- Electrical Requirements: 100/120V 50/60Hz
- Weight (English): 214 lb.
- Depth (English) Interior: 27.5 in.
- Vial Capacity: 5500 to 6318

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com