

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Warning: file_put_contents(/home/beqsaleshost/public_html/logs/deprecated.php): Failed to open stream: Permission denied in /home/beqsaleshost/public_html/libraries/vendor/joomla/filesystem/src/File.php on line 290

Thermo Scientific Double Glass Door Lab Refrigerator



Capacity: 51.1 cu. ft. (1447 L)

Rating: Not Rated Yet

Price

Sales price \$5,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Thermo Scientific TSX5005CA is a high-performance chromatography refrigerator with double glass doors, designed for precise temperature control (2°C to 8°C) and easy access to instruments in labs, healthcare, and research settings. It features V-drive technology for energy efficiency, temperature uniformity, and reliable storage of samples, reagents, and chromatography equipment.

Key Features

Capacity: 51.1 cu. ft. (1447 L) with 4 full + 2 half epoxy-coated adjustable shelves for flexible storage.

Temperature Range: 2°C to 8°C (setpoint ~5°C) with excellent uniformity and rapid recovery, plus >450W reserve capacity for heat loads from instruments.

Double Glass Doors: Self-closing with 90° stop, allowing easy viewing and access while minimizing temperature fluctuations; includes ergonomic handles and standard door locks.

V-Drive Technology & Sustainability: Variable-speed compressor adapts to usage patterns for significant energy savings; uses natural R290

Refrigerators: Thermo Scientific Double Glass Door Lab Refrigerator

refrigerant; low noise (~52 dBA); ACT-certified and GMP/ISO cleanroom compatible.

User-Friendly & Safety Features: Microprocessor control with digital display, alarms, setpoint security, bright LED interior lighting, two 2" access ports, hydrocarbon safety switch-protected internal duplex outlet, and connectivity options (e.g., DeviceLink, remote alarms).

Dimensions (Exterior): ~37.8" D x 56.5" W x 78.5" H; operates on 115V/60Hz.

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.

226 Miller Ave.

South San Francisco, CA 94080

650-871-5707 office

650-276-7487 fax

sales@biotechequipmentsales.com