


Thermo Scientific EC1000XL Power Supply



Current: 500 mA

Rating: Not Rated Yet

Price

Sales price \$850.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Ideal for DNA, RNA and protein electrophoresis, as well as electroblotting. Simple to use, multi-purpose, stackable, and lightweight. Power supply includes soft-touch key pad and four jacks.

Specifications

Current: 500 mA

Model EC1000XL

Max. Current: 500 mA

Timer: 0 to 99 hr. 59 min.

Display: Voltage, current or power

No. of Jacks: 4

Dimensions (L x W x H): 12 x 10.75 x 5.25 in. (31 x 27 x 13.3 cm)

Electrical Requirements: 120 V 50/60 Hz

Electrophoresis: Thermo Scientific EC1000XL Power Supply

Max. Voltage: 1000 V

Comes with our 90-day warranty.

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

Biotech Equipment Sales, Inc.
226 Miller Ave.
South San Francisco, CA 94080
650-871-5707 office
650-276-7487 fax
sales@biotechequipmentsales.com