


## Thermo Scientific Fiberlite F9-6x1000 LEX Fixed Angle Rotor



9000 rpm; 17,600 x g for LYNX centrifuge

Rating: Not Rated Yet

**Price**

Sales price \$11,500.00

[Ask a question about this product](#)

Manufacturer [Thermo Scientific](#)

Description

Fiberlite F9-6x1000 LEX Fixed Angle Rotor for Use With LYNX 6000 Superspeed Centrifuge

Specifications

Angle: 20°

K Factor: 3,415

Max. Tube Diameter: 98 mm

Max. Tube Length: 195 mm

Net Weight: 34.5 lb. (15.7 kg)

Radius (Metric): 194 mm

Radius (Metric) Min.: 65 mm

Biocontainment: Yes

Capacity: 6 x 1000 mL

Material: Carbon Fiber Composite

Max. RCF: 17,568 x g

Max. Speed: 9,000 rpm

Type: Fixed Angle Rotor

*Comes with our 90-day warranty.*

Biotech Equipment Sales, Inc.: based in South San Francisco, the Biotech capital of the world, since 2009. We are dedicated to bringing high-quality pre-owned laboratory instrumentation to our scientific community.

**Biotech Equipment Sales, Inc.**  
226 Miller Ave.  
South San Francisco, CA 94080  
650-871-5707 office  
650-276-7487 fax  
[sales@biotechequipmentsales.com](mailto:sales@biotechequipmentsales.com)